




Fronius gas cooled robotic PAP Hose pack with open connector

Example: Fronius gas cooled robotic PAP




Hose pack with a miller power pin









Fronius gas cooled robotic PAP Hose pack mounted on a miller wire feeder



Powerpins part numbers


PowerPin	Article-Number	Suitable for ...
	42,1000,0000	Lincoln
	42,1000,0001	Miller
	42,1000,0002	Hobart

Installing PowerPin to MHP Robot Hosepacks



Clamping piece

- You have to use the same like for Euro.
- Tool for mounting is in the power pin included
- 44,0350,3703 Clamping piece 0,8/E
- 44,0350,3704 Clamping piece 1,0/E
- 44,0350,3705 Clamping piece 1,2/E
- 44,0350,3706 Clamping piece 1,6/E
- We have a second version too, without thread for fixing, that you can change the inner liner from the contact tip side (FL means Front Loaded)
- 44,0350,0005 Clamping piece PAP 1,0 FL Euro
- 44,0350,0008 Clamping piece PAP 1,2 FL Euro
- 44,0350,0039 Clamping piece PAP 1,6 FL Euro





All information is without guarantee in spite of careful editing – liability excluded.

Intellectual property and copyright: all rights reserved.
Copyright law and other laws protecting intellectual property apply to the content of this presentation and the documentation enclosed (including texts, pictures, graphics, animations etc.) unless expressly indicated otherwise. It is not permitted to use, copy or alter the content of this presentation for private or commercial purposes without explicit consent of Fronius.